

Standard Operating Procedure for the
Analysis of Silica in Fresh Waters
CCAL 32A.0

Cooperative Chemical Analytical Laboratory
Forestry Sciences Laboratory
Oregon State University
3200 SW Jefferson Way
Corvallis, Oregon

Prepared by Kathryn Motter and Cameron Jones
March 2006

**Standard Operating Procedure for the
Analysis of Silica In Fresh Waters
CCAL 32A.0**

Table of Contents

1.0	Scope and Application	3
2.0	Summary of Method	3
3.0	Definitions	3
4.0	Interferences	4
5.0	Safety	4
6.0	Equipment and Supplies	4
7.0	Reagents and Standards	5
7.1	Preparation of Reagents	5
7.2	Preparation of Standards	6
8.0	Sample Handling and Storage	6
9.0	Quality Control	6
10.0	Calibration and Standardization	6
11.0	Procedure	7
11.1	Calibration and Analysis Procedure	7
11.2	System Notes	7
12.0	Data Analysis and Calculations	8
13.0	Method Performance	8
14.0	Pollution Prevention	8
15.0	Waste Management	8
16.0	References	9
17.0	Tables, Diagrams, Flowcharts, and Validation Data	9
17.1	Ammonia Nitrogen Reaction Manifold	9
17.2	Data System Parameters	10

Standard Operating Procedure for the Analysis of Silica In Fresh Waters CCAL 32A.0

1.0 Scope and Application

- 1.1 This method details the determination of silicate in fresh waters by automated colorimetric analysis. The practical range of determination for this method is 0 to 20 mg SiO₂/L. Method detection limit for this analysis is 0.1 mg/L as silica.

2.0 Summary of Method

- 2.1 An automated analysis method is used for the colorimetric determination of silica in fresh water. Silica reacts with molybdate in an acidic environment to form yellow beta-molybdosilicic acid. The acid is reduced with ascorbic acid to form a heteropoly molybdous-blue complex that absorbs at 660 nm (See Reference 16.4). Concentration of silica is determined by comparison of absorbance signal with calibration results obtained from prepared standards of varying concentrations.

3.0 Definitions

- 3.1 DI water: Water that has been through a deionization system to produce water similar to ASTM Type I reagent with 16.7 Mohms resistivity (ASTM) (Reference 16.3).
- 3.2 Method Detection Limit (MDL): The minimum concentration of an analyte that can be measured and reported with 99% confidence, based on a one-sided 99% confidence interval (t -value at a significance level of 0.01 and $n-1$ degrees of freedom) from at least seven repeated measurements of a low concentration standard measured within an analysis run.

$$\text{MDL} = ts$$

Where,

t = Student's t value at a significance level of 0.01 and $n-1$ degrees of freedom

s = standard deviation of at least seven repeated measurements of a low level standard

4.0 Interferences

- 4.1 Oxalic acid is added to the sample stream before addition of ascorbic acid to eliminate interference from phosphates.
- 4.2 Tannin and large amounts of iron, color, turbidity and sulfide interfere.
- 4.3 Freezing samples may reduce silica concentration (Reference 16.1)
- 4.4 Samples, standards and reagents should be stored in polyethylene containers.
- 4.5 Turbidity and color that persists in the sample after filtration may absorb at the 660 nm. Reanalyze samples without the color reagent as turbidity/color blanks. Any response from the turbidity/color blanks must then be subtracted from the initial analysis response.

5.0 Safety

- 5.1 The toxicity or carcinogenicity of each reagent has not been precisely determined; however, each chemical should be regarded as a potential health hazard. Exposure to these chemicals should be reduced to the lowest possible level. Cautions are included for known extremely hazardous materials.
- 5.2 The following chemicals have the potential to be highly toxic or hazardous. For detailed explanations, consult the MSDS.
 - 5.2.1 Sulfuric acid
 - 5.2.2 Oxalic acid

6.0 Equipment and Supplies

Note: *Brand names, suppliers and part numbers are for illustrative purposes only. No endorsement is implied. Equivalent performance may be achieved using apparatus and materials other than those specified here, but demonstration of equivalent performance that meets the requirements of this method is the responsibility of the laboratory.*

- 6.1 Balance measuring to at least 0.1 mg
- 6.2 Volumetric flasks and pipettes as required

- 6.3 Technicon Auto-Analyzer II
 - 6.3.1 Multichannel proportioning pump
 - 6.3.2 Colorimetric Detector
 - 6.3.3 Data system
 - 6.3.4 Alpkem manifold and method 105-71W/B
- 6.4 Safety glasses
- 6.5 Nitrile gloves
- 6.6 Lab coat or apron
- 6.7 Laboratory Exhaust Fume Hood
- 6.8 High Density Polyethylene (HDPE) bottles

7.0 Reagents and Standards

7.1 Preparation of Reagents

Solutions 7.1.2 through 7.1.3 need to be filtered through prewashed GF/F filters after preparation

- 7.1.1 *0.1 N Sulfuric Acid:*
Carefully add 2.8 mL concentrated sulfuric acid to 800 mL DI water. Cool and dilute to 1 L.
- 7.1.2 *Ammonium Molybdate:*
Dissolve 10 g of ammonium molybdate in 1 L of 0.1 N sulfuric acid. Store in amber polyethylene container. Solution is stable for 2 weeks.
- 7.1.3 *Oxalic Acid:*
Dissolve 50 g of oxalic acid in 900 mL of DI water. When dissolution complete, dilute to 1 L. Store in amber polyethylene container. Reagent is stable.
- 7.1.4 *15 % Sodium Laurel Sulfate:*
Dissolve 4.41 g sodium laurel sulfate in 25 mL DI water.
- 7.1.5 *Ascorbic Acid:*
Dissolve 17.6 g pf U.S.P. quality ascorbic acid in 500 mL of DI water containing 50 mL acetone. Add 20 mL of 15 % sodium laurel sulfate, and dilute to 1 L with DI water.

7.2 Preparation of Standards

7.2.1 A 1000 mg/L commercial silica stock standard is used. Working standards in concentrations of 25.0, 15.0 and 5.0 mg/L silica are prepared from the stock standard.

8.0 Sample Handling and Storage

8.1 Unfiltered samples are filtered upon receipt through glass fiber filters into clean HDPE bottles and stored at 4°C in the dark. Samples are analyzed within 48 hours to ensure sample integrity. If samples must be held prior to analysis, they are stored refrigerated at 4°C.

9.0 Quality Control

- 9.1 Preparation of standards is recorded on worksheets and documented by weight of standard added to a given flask before dilution to volume with DI water. All records of certification are kept on file at CCAL Laboratory.
- 9.2 Blank: DI water run after the calibration and before and after each check standard.
- 9.3 Quality Control Check Standard: Calibration standards run in rotation every ten samples to monitor stability and validate the calibration.
- 9.4 Method Detection Limit (MDL): Established for each analyte. Based on a one-sided 99% confidence interval (t-value) from at least seven repeated measurements of a low concentration standard. The t-distribution value is multiplied by the standard deviation of the population (n-1) to obtain the MDL.
- 9.5 Analytical Duplicate: Separate analysis from the same sample aliquot. Run a minimum of once every ten samples.

10.0 Calibration and Standardization

- 10.1 Balances: calibrated yearly by external vendor.
- 10.2 Pipette delivery checked by weight to within 2% of theoretical weight of aliquot volume.

- 10.3 Calibration curve with $r^2 \geq 0.998$. (See 17.2 for calibration data set-up.)
- 10.4 Calibration verification with check standards, monitored throughout the run. If measurement exceeds +/- 10% of the theoretical value, the analysis should be terminated and the instrument recalibrated. The calibration must be verified before continuing analysis.

11.0 Procedure

11.1 Calibration and Analysis Procedure

- 11.1.1 Prepare reagents and standards as outlined in Section 7.
- 11.1.2 Set up manifolds as shown in Section 17.1.
- 11.1.3 Samples are injected into the reaction path at a fixed time interval, determined by the cam timing set in the auto-sampler. Setup or confirm data system parameters as detailed in Section 17.2.
- 11.1.4 Pump DI water through all reagent and sample lines. Check for leaks and stable bubble pattern (smooth flow). Pump reagents through all lines until the system equilibrates (minimum of one hour).
- 11.1.5 Record sample id's in the data template.
- 11.1.6 Calibrate the instrument with standards. Calibration regression equations must have $r^2 \geq 0.998$.

11.2 System Notes

- 11.2.1 If the baseline is excessively noisy, clean the manifold using the following procedure:
- Place all reagent and carrier lines in rinse water and pump to clear reagents.
 - Place all lines in 1.2N HCl solution for several minutes.
 - Place lines in DI water and pump until thoroughly rinsed.
- 11.2.2 If baseline has a large number of air spikes, check pump tubes for excessive wear and replace as necessary.
- 11.2.3 Poor sensitivity or insufficient color development may be the result of old ammonium molybdate solution. Prepare fresh and flush system thoroughly.
- 11.2.4 Turbidity and color that persists in the sample after filtration may absorb at the 660 nm. Reanalyze samples without the color reagent as turbidity/color blanks. Any response from the turbidity/color blanks must then be subtracted from the initial analysis response.

12.0 Data Analysis and Calculations

- 12.1 The data system prepares a calibration curve by plotting response of injected standards versus known standard concentration. The resulting regression equation is used to calculate the sample concentration.
- 12.2 All results and sample information are filed in the analysis data system by analysis run. Details specific to the instrumental analysis are noted in the Instrument Run Log created and maintained for the AAIL. Analytical results are entered into electronic format and entries are verified by a second person.
- 12.3 Historically, CCAL has reported results as mg Si/L. To convert from SiO₂ to Si, divide results by 2.137; conversely, to convert from Si to SiO₂, multiply by 2.137.

13.0 Method Performance

- 13.1 This method was validated through inter-laboratory studies. The CCAL participates in the USGS Standard Reference Water QA program.

14.0 Pollution Prevention

- 14.1 The chemicals used in this method pose little threat to the environment when properly managed.
- 14.2 All standards and reagents should be prepared in volumes consistent with laboratory use to minimize the volume of disposable waste.
- 14.3 For further information on pollution prevention consult *Less is better: Laboratory Chemical Management for Waste Reduction*, available from the American Chemical Society's Department of Government Relations and Science Policy, 1155 16th Street NW, Washington D.C. 20036, (202) 872-4477.

15.0 Waste Management

- 15.1 It is the laboratory's responsibility to comply with all federal, state and local regulations governing waste management, and to protect the environment by minimizing and controlling all releases from fume hoods

and bench operations. Compliance with all sewage discharge permits and regulations is required.


- 15.2 For further information on waste management, consult “The Waste Management Manual for Laboratory Personnel”, and “Less is Better: Laboratory Chemical Management for Waste Reduction”, both available from the American Chemical Society’s Department of Government Relations and Science Policy, 1155 16th Street NW, Washington DC, 20036.

16.0 References

- 16.1 Standard Methods For The Examination of Water and Wastewater, Method 4500-SiO₂, Automated Method for Molybdate-Reactive Silica. American Public Health Association. 21st Edition, 2005.
- 16.2 Code of Federal Regulations. Protection of Environment. Section 40, Appendix B to Part 136. Definition and procedure for the determination of the method detection limit. Revision 1.11. Revised July 1, 1990. Office of the Federal Register, National Archives and Records.
- 16.3 ASTM. American Society for Testing and Materials. Standard Specifications for Reagent Water. D1193-77 (Reapproved 1983). Annual Book of ASTM Standards, Vol. 11.01. ASTM: Philadelphia, PA, 1991.
- 16.4 Technicon AutoAnalyzer II; Industrial Method No. 105-71W/B; rev Jan 1976. Technicon Industrial Systems; Tarrytown, N.Y. 10591

17.0 Tables, Diagrams, Flowcharts, and Validation Data

- 17.1 Silica Reaction Manifold



Silica Manifold Specifications

Carrier is molybdate reagent

Interference Filter is 660 nm

Pump tubing is Tygon

Manifold tubing is 0.030 mm i.d.

Mixing coils are 20 turn with mid injection fitting, 20 turn with front injection, and 22 turn.

15 mm x 2 mm flow cell

17.2 Data System Parameters

Cycle throughput: 30 samples/hr

Cycle Period: 120 s

Analyte Data:

Concentration Units mg SiO₂/L

Calibration Data:

Level	1	2	3
Concentration mg/L	25.0	15.0	5.0

Calibration Fit Type: Linear Regression

Sampler Timing:

Minimum Probe in Wash Period: 30 s

Probe in Sample Period: 90 s